

OPHTHALMOLOGY INSTRUMENTS CATALOG



Dental | Beauty | Veterinary



 **VOLGO CARE**

Barraquer

Wire Speculum

K1-5000 standard wire, 10mm blades

K1-5010 standard wire, 15mm blades

K1-5009 heavy wire, 15mm blades




Grandon-Barraquer

Wire Speculum

with suture notches

Used to obtain additional traction for the upper eyelid.

K1-5011 15mm blades



Temporal Approach


Wire Speculum

with spring on nasal side

K1-5012 15mm closed blades

K1-5015 15mm open blades

K1-5014 15mm solid blades




Barraquer

Wire Speculum

solid blades

K1-5020 10mm blades

K1-5030 15mm blades



Speculums



Feaster

Wire Speculum

K1-5032 20mm solid blades



Kratz-Barraquer


Wire Speculum

15mm blades

one solid and one open wire

K1-5044 for right eye

K1-5045 for left eye



Kratz-Barraquer

Wire Speculum


with open blades

For unobstructed access to the cataract incision with phaco or I/A tip.

K1-5049 heavy wire, 15mm blades

K1-5050 standard wire, 15mm blades

K1-5051 standard wire, 13mm blades




Kratz-Barraquer

Wire Speculum
with wings

*For retracting the lids
up and away from the globe
to eliminate external pressure.*

K1-5060 15mm blades





McIntyre

Wire Speculum
gently curved springs

*V-shaped blades are designed to
retain surgical drape under the eyelids.*

K1-5070 closed wire, 13mm blades

K1-5071 open wire, 13mm blades




Nevyas

Wire Speculum

*Designed to retain surgical drape under the
eyelids while allowing unobstructed access to
the cataract incision with phaco or I/A tip.*

K1-5073 16mm blades



Speculums



Lieberman

Aspirating Speculum
adjustable mechanism
Kratz style 15mm blades with aspirating ports to remove excess fluid from the surgical field.
Supplied with silicone tubing and luer-lock adaptor.

K1-5171 complete

K1-5179 replacement silicone tubing, 100cm (39")



Lieberman

Aspirating Speculum
adjustable mechanism
V-shaped 15mm blades with aspirating ports to remove excess fluid and retain the surgical drape under the lid.
Supplied with silicone tubing and luer-lock adaptor.

K1-5175 complete

K1-5179 replacement silicone tubing, 100cm (39")



Sauer

Lid Speculum

K1-5300 for infants, 11mm blades

K1-5302 for premature infants, 2mm blades




Frankel

Infant Lid Speculum
adjustable
*For precise control of
infant eyelid retraction.*
K1-5310 10mm blades



Alphonso

Infant Lid Speculum
finger grips and
interlacing blades to
facilitate insertion
K1-5340




Infant


Lid Speculum
closed wire blades
gently curved springs
K1-5350 4mm blades



Tanna-Lieberman
Eye Speculum
adjustable mechanism
wire blades with traction
suture anchor posts
K1-5668 15mm blades



Lieberman
Eye Speculum
Kratz style open wire blades
adjustable mechanism
K1-5670 10mm blades
K1-5671 15mm blades



Lieberman
Eye Speculum
adjustable mechanism
*V-shaped blades are designed
to retain the surgical drape
under the eyelids.*
K1-5675 15mm blades



Scleral Shield

stainless steel

non-reflective exterior surface

To protect the eye during blepharoplasty and skin resurfacing procedures.

K1-8480 1 small, 1 large (set)

K1-8481 20 x 28mm, small (pair)

K1-8482 22 x 30mm, large (pair)



Jaeger

Lid Plate

clear plastic

K1-8500



Jaeger

Lid Plate

stainless steel

K1-8520



Wright

Fascia Needle

1mm x 6mm oval hole

140mm long

K1-8540



Retractors



Knapp
Lacrimal Retractor
4 prongs
K1-9000 blunt
K1-9001 sharp



Walton
Conjunctiva Retractor
6mm wide
thin solid blade
K1-9008



Helveston
Tissue Retractor
thin curved blade
K1-9010 "Barbie" 7mm wide
K1-9012 "Big Barbie" 9mm wide
K1-9019 "Great Big Barbie" 11mm wide concave blade



K1-9019

Tooke
Corneal Knife
straight blade
with curved cutting edge
K2-3650



Corneal Dissector
angled, straight blade
For intrastromal dissection.
K2-3660



Corneal Dissector
angled, curved blade
For intrastromal dissection.
K2-3661



Gill
Corneal Knife
curved cutting edge
K2-3700 straight
K2-3710 light curve



Knives • Spuds



Fukasaku

"Hockey" Knife
sharp 3mm front edge
semi-sharp 6mm bottom edge
K2-3735




Golf Club

Foreign Body Spud
angled, lightly rounded tip
K2-4000




Francis

Foreign Body Spud
angled, flag shaped tip
K2-4050



Ellis


Foreign Body Spud
lightly curved and rounded tip
K2-4100





Meyerhoefer
Chalazion Curette

- K3-1000** size 0, 1.5mm dia.
- K3-1010** size 1, 1.75mm dia.
- K3-1020** size 2, 2.25mm dia.
- K3-1030** size 3, 2.5mm dia.
- K3-1040** size 4, 3.5mm dia.



Nightingale
Posterior Capsule
Polishing Curette
1.75mm diameter open ring
with semi-sharp edge
K3-1100



Shepherd
Capsule Polishing Curette
1.75mm diameter low profile
open ring semi-sharp edges
on top and bottom
*For polishing posterior and
anterior capsule.*
K3-1102



Quality without Compromise


Curettes & Spatulas



De La O
Nucleus Delaminator
1.75mm diameter ring
K3-1108



Blumenthal
Conjunctiva Dissector
with blunt disc shaped tip
K3-1124



Barraquer
Corneoscleral Brush
sable, with plastic handle
K3-1750



Buratto
LASIK Flap Protector
*For retracting and
protecting nasally hinged
corneal flaps.*
K3-1770



Eye Instruments

Buratto

LASIK Flap Protector
 blades angled 15°
*For retracting and protecting
 superiorly hinged corneal flaps.*
K3-1772



Mannis-Buratto

LASIK Flap Protector
 blades angled 45°
*For retracting and protecting
 superiorly hinged corneal flaps.*
K3-1773



Rowen

Combination LASIK Flap
 Protector and Spatula
K3-1774



Camellin Style

LASEK Alcohol Well
*For retaining alcohol solution
 on the corneal epithelium.*
K3-1820 8.5mm diameter
K3-1822 9.5mm diameter

BOTTOM VIEW



Quality without Compromise

Spatulas



Sloane

LASEK Micro Hoe
semi-circular tip with
beveled edge
*For lifting the precut margin
of the epithelial flap.*

K3-1845



Sloane

LASEK Epi Peeler
semi-sharp edges
*For peeling and folding back
the epithelial flap.*

K3-1845



Sloane


LASEK Micro Hoe &
Epi Peeler
K3-1848



Sloane

LASEK Epi-Dissector
spatulated tip with
semi-sharp edges
For undermining the epithelium.

K3-1850



Sloane

LASEK Flap Repositor
 gently curved spatulated tip
 with beveled smooth edges
*For repositing and smoothing
 the epithelial flap after laser
 ablation.*

K3-1855



Culler

Iris Spatula
K3-2100 1mm wide
K3-2110 2mm wide



Wheeler

Double-End Spatula
 lightly curved
K3-2120



Barraquer

Iris Spatula
 round, smooth
K3-2300 0.5mm diameter
K3-2310 0.25mm diameter



Spatulas



Gimbel

Nucleus Spatula
round, smooth
angled 30° at
10mm from the tip
0.25mm diameter
K3-2312



Escaf

Nucleus Manipulator
angled 45° at
13mm from the tip
K3-2347



Knolle

Nucleus Spatula
malleable stainless steel, flat
0.5mm wide
K3-2350



Koch

Nucleus Spatula
gently curved duckbill
shaped tip
K3-2354 smooth
K3-2355 with notches



Nucleus Manipulators

Akahoshi

Nucleus Sustainer
 Designed to provide
 counterpressure when
 used with the
 Akahoshi PreChopper.
K3-2357



Escaf

Nucleus Sustainer
 smooth ball tip
 textured inferior surface
K3-2359



Lieberman

MicroFinger
 Nucleus Manipulator
K3-2360 for use in right hand
K3-2361 for use in left hand



Lieberman

MicroFinger and Spatula
K3-2362 for use in right hand
K3-2363 for use in left hand




Quality without Compromise

Nucleus Choppers



Fukasaku-Lieberman

Phaco Spatula
a combination of the
small pupil snapper hook
and the MicroFinger
K3-2364 for use in right hand
K3-2365 for use in left hand



Masket

Phaco Spatula
a combination of the Rosen
Phaco Splitter and a 12mm
long round spatula with disc
shaped tip
K3-2367



Chang

Phaco Chopper
a combination of a modified
MicroFinger and a 0.9mm
long "Quick Chop" tip
K3-2368 for use in right hand
K3-2369 for use in left hand



Minardi

Phaco Chopper
0.7mm long
"Quick Chop" tip
K3-2370



Langerman

Phaco Chopper
elongated Kuglen style
push-pull tip with
wedge shaped superior
and inferior edges
For bi-directional chopping.
K3-2372



Hwang

Phaco Chopper
1mm long
vertical chopper
and a chisel
tipped spatula
K3-2374



Agarwal

Karate Phaco Chopper
semi-sharp, wedge shaped
inferior edge with a
pointed tip
K3-2375



Quality without Compromise

Nucleus Choppers



Chang

MicroFinger Chopper
curved finger shaped tip
with defined inner edge
for horizontal chopping
K3-2377



Warren

Canine Phaco Chopper
a combination of a 3.5mm
long MicroFinger and a 2mm
long pointed "Quick Chop" tip
For veterinary use only.
K3-2379



Brierley

Nucleus Splitter
with incision sealing shaft
K3-2380



Dodick-Kammann

Nucleus Chopper
1.5mm long tip

For bimanual nucleus division.

K3-2382 straight

K3-2383 angled at 12mm



Dodick-Kammann

Nucleus Chopper
2mm long tip

For bimanual nucleus division.

K3-2384 straight

K3-2385 angled at 12mm



Escaf

Nucleus Divider
1.25mm long

wedge shaped tip

For vertical chopping.

K3-2387



Nucleus Choppers



Neuhann
Nucleus Divider
strongly curved spatulated tip
K3-2388



Nagahara
Nucleus Spatula
flat inferior edge
K3-2390



Fukasaku
Snapper Hook
wedge shaped inferior edge
*For the "Snap & Split"
phaco technique.*
K3-2392



Fukasaku
Small Pupil
Snapper Hook
wedge shaped inferior
edge with blunt tip and
front notch for engaging iris
K3-2393



Eye Instruments

Rosen
Phaco Splitter
wedge shaped inferior edge
with blunt tip
K3-2395



Rosen
Phaco Splitter
60° angled wedge shaped
inferior edge with blunt tip
K3-2396



Rowen
Phaco Spatula
a combination of the
Rosen Phaco Splitter
and the Bechert
Nucleus Rotator
K3-2398



Bechert
Nucleus Rotator
blunt forked tip
angled at 7mm
K3-2400



Nucleus Manipulators

Jaffe-Bechert
Nucleus Rotator
blunt forked tip angled at 9mm
K3-2410



Clayman
Nucleus Rotator
0.5mm wide flat spatula
with V-notched tip
K3-2415



Drysdale
Nucleus Manipulator
paddle shaped tip
K3-2418




Drysdale
Nucleus Manipulator
and Anterior Capsule Polisher
paddle shaped tip
with defined edges
on anterior surface
K3-2419



Eye Instruments

VOLGO CARE Nucleus Manipulators

Pisacano
Nucleus Rotator
horizontal flat disc
shaped tip
K3-2420




Pisacano
Nucleus Rotator
vertical flat disc
shaped tip
K3-2421



Siepser
Endocapsular Controller
trumpet shaped tip
K3-2425



Anis
Ball Tip
Nucleus Manipulator
with 0.5mm diameter
ball tip
K3-2426 straight
K3-2427 angled



Quality without Compromise

Manipulators & Spatulas



Tennant
Nucleus Rotator
with 0.6mm diameter
ball shaped tip
K3-2428



Lippman
Nucleus Spatula
0.5mm diameter round
shaft with blunt forked tip
K3-2430



Girard
Synechia Spatula
0.5mm wide with
semi-sharp tip
K3-2500



Obstbaum
Synechia Spatula
0.25mm wide with
semi-sharp tip
K3-2520




Moran

LASIK Enhancement
Spatula
Fukasaku spatula
on one end and
miniature Paton
spatula on the
opposite end
K3-2524



Hersh

IntraLase® Flap Spatula
beveled tips
a combination of a
6mm spatula and
a 10mm spatula
K3-2526



Binder

IntraLase® Flap Elevator
1mm wide
12mm long
semi-sharp edges
K3-2528



Fukasaku

LASIK Spatula
rounded anterior and
flat posterior surface
with beveled tip
*For lifting and dissecting
the corneal flap.*
K3-2530



LASIK Spatulas




Guimaraes

LASIK Spatula
gently curved, bi-convex
semi-sharp edges
K3-2532



Katzen

LASIK Flap "Unzipper"
combines a modified
Sinsky Hook and
a bi-convex spatula with
beveled notches
For LASIK enhancements.
K3-2534



Maddox

LASIK Spatula
combines a flat spatula and
a cylindrical spatula
*For elevating and
repositing the flap.*
K3-2535



LASIK

Flap Manipulator
delicate olive tip
*Designed to locate
and open a previously
made LASIK flap.*
K3-2540



Eye Instruments

Hirschman
 Lens Spatula
 flat 0.5mm wide
 with multiple notches
K3-2550



Jaffe
 Lens Spatula
 flat, 0.5mm wide
K3-2600



Lester
 IOL Manipulator
 0.25mm diameter tip
Ideal for multi-directional manipulation.
K3-2680 straight, with guard
K3-2690 angled



IOL Manipulators



Lester

IOL Manipulator
0.25mm diameter tip
vaulted shank
*For engaging the inferior haptic
without touching the dome
of the IOL.*
K3-2691



Maloney

IOL Manipulator
0.25mm diameter tip
vaulted shank
with guard
K3-2692



Hunkeler

IOL Manipulator
with 0.3mm diameter
hemi-ball tip
K3-2695



Alio-Kelman

IOL Hook
0.2mm diameter tip
K3-2697



VOLGO CARE Manipulators & Spatulas

Aker

Lens Pusher
smooth, fish mouth
shaped tip
with guard
K3-2720



Bechert

Lens Pusher
vertical Y-shaped tip
angled
K3-2727




Elschnig

Cyclodialysis Spatula
flat, 1mm wide
K3-3000



Castroviejo

Synechiaie Spatula
0.5mm wide with
semi-sharp tips
K3-3150




Quality without Compromise

Spatulas & Loops




Castroviejo
Cyclodialysis Spatula
K3-3200 0.25mm wide
K3-3210 0.5mm wide



Kimura
Platinum Spatula
K3-3500



Keener
Lens Divider
with three preformed
snare wire loops
K3-3900 complete
K3-3910 wire loops (pk/3)



Lewis
Lens Loop
K3-4000 small
K3-4010 large



New Orleans
Lens Loop
K3-4030



Wilder
Lens Loop
serrated
K3-4100



Kirby
Lens Loop and
Muscle Hook
K3-4105



Lens Loops



Morrison
Lens Loop and Probe
K3-4120



Kansas
Lens Loop
serrated, 3mm wide
*For manual small
incision nucleus
fragment removal.*
K3-4140



Keener-Artt
Lens Loop
serrated, 3mm wide
strong curve
*For manual small
incision nucleus
fragment removal.*
K3-4141



Alfonso
Nucleus Trisector
1.5mm wide
Used with K3-4156
Alfonso Cutting Platform
For manual small incision surgery.
K3-4145



VOLGO CARE Nucleus Fragmentation

Bucher-Kansas

Nucleus Trisector
2.5mm wide
Used with K3-4157
Alfonso Cutting Platform
For manual small incision surgery.
K3-4146



Kansas

Nucleus Trisector
4mm wide
Used with K3-4152
Kansas Vectis
For manual small incision surgery.
K3-4148



Kansas

Nucleus Bisector
For manual small incision surgery.
K3-4150



Kansas

Nucleus Vectis
Used as a cutting block for
nucleus bisector and trisector.
For manual small incision surgery.
K3-4152



Quality without Compromise

Nucleus Fragmentation

Alfonso

Cutting Platform

Used as a cutting block for nucleus trisector.

For manual small incision surgery.

K3-4156 2mm wide, used with K3-4145

K3-4157 3mm wide, used with K3-4146



Giannetti

Nucleus Capture Spoon

K3-4160



Giannetti

Nucleus Capture Spatula


K3-4162



Giannetti

Nucleus Capture Hook

K3-4164



Eye Instruments

Culler

Lens Spoon
sterling silver
K3-4200



Paton

Spatula & Spoon
For corneal button
transfer.
K3-4255



Rosenwasser

Lamella Donor Shovel
with retaining rim on 2 sides
For atraumatic insertion
of the donor lamella.
K3-4260 4mm wide



Rosenwasser

Lamella Donor Shovel
with retaining rim on 3 sides
For atraumatic insertion
of the donor lamella.
K3-4263 5mm wide



Spoons




Bunge
Evisceration Spoon
K3-4400 small
K3-4410 large



Wells
Enucleation Spoon
K3-4450



Graether
Collar Button
micro iris retractor
K3-4900 straight,
with guard
K3-4910 angled



Eye Instruments

Graether

Irrigating Collar Button
micro iris retractor and
posterior capsule polisher
K3-4920 straight
K3-4930 angled



Keuch

Small Pupil Dilator
with guard
*For enlarging a small pupil
through the sideport and
primary incision.*
K3-4950



Iris Retractors

six 10mm long flexible hooks
in an autoclavable case
reusable
K3-4970

*Also available
sterile/disposable
see page 204*



Hooks



Z-Hook

Manipulator
Z-shaped blunt tip
0.15mm diameter
*For manipulation
of the donor lamella.*

K3-5000



Reversed

Sinskey Hook
0.2mm diameter blunt tip
*For scoring the recipient
endothelium and positioning the
donor lamella.*

K3-5002



Bonn

Micro Iris Hook
extra delicate
blunt, round handle
with guard

K3-5050



Bonn

Micro Iris Hook
blunt, flat handle

K3-5060



Shepard
Iris Hook
blunt
K3-5110



Jaffe
Lens Manipulating Hook
0.15mm diameter blunt tip
straight
K3-5160



Clayman-Jaffe
Lens Manipulating Hook
0.15mm diameter blunt tip
angled
K3-5170



Sinsky
Lens Manipulating Hook
0.25mm diameter blunt tip
straight
K3-5180



Quality without Compromise

IOL Manipulators




Harris-Sinskey

Lens Manipulating Hook
0.25mm diameter blunt tip
angled
K3-5190




Sinskey

Lens Manipulating Hook
round handle
0.25mm diameter blunt tip
K3-5192 straight, with guard
K3-5193 angled



Nevyas

Soft IOL Manipulator
0.25mm diameter
textured tip
straight, with guard
K3-5194



Sinskey II

Lens Manipulating Hook
0.2mm diameter blunt tip

K3-5220 straight

K3-5230 angled



Sinskey II

Lens Manipulating Hook
0.2mm diameter blunt tip
round handle

K3-5232 straight, with guard

K3-5233 angled



Lewicky

Lens Manipulating Hook
0.15mm diameter blunt tip
angled

K3-5250



Hirschman

Iris Hook
blunt, angled

K3-5400



Quality without Compromise

Kuglen Hooks



Kuglen
Iris Hook and Lens
Manipulator
push-pull model
K3-5500 straight
K3-5520 angled



Kuglen
Iris Hook and Lens
Manipulator
clover leaf model
K3-5510



Kuglen
Iris Hook and Lens
Manipulator
push-pull model
round handle
K3-5522 straight, with guard
K3-5523 angled



Muscle Hooks



Hofmann-Thornton
Globe Fixation Ring
with swivel handle
16mm diameter ring
multiple blunt teeth
K3-6168



Stevens
Tenotomy Hook
blunt, small
K3-6600



Helveston
"Teaser" Hook
extra delicate, blunt
K3-6610



Graefe
Muscle Hook
K3-6710 size 1
K3-6720 size 2
K3-6730 size 3



Gass
Retinal Detachment Hook
with oval hole
K3-6740



Helveston
Finder Hook
extra delicate with lightly
angled cone shaped tip
K3-6760 large
K3-6761 small



Helveston
Muscle Hook
extra delicate, with blunt
spatulated tip
K3-6778 8mm
K3-6780 10mm
K3-6782 12mm



Muscle Hooks



Jameson
Muscle Hook
K3-6800 small
K3-6810 large



Du Plessis
Muscle Hook
double arm
K3-6814



Guyton
Small Incision
Muscle Hook
K3-6820



Green
Muscle Hook
K3-6850



Calipers



Castroviejo
Caliper
20mm spread
K3-9000 straight
K3-9002 angled
K3-9010 curved



Osher
Incision Caliper
*For internal measurement
of cataract incisions from
1 to 10mm.*
K3-9012



Helveston
Scleral Ruler
15mm long, with notches
in 5mm increments
single marking tooth at tip
*For measuring and marking
the distance from the limbus
to the location of the muscle
reinsertion.*
K3-9030



Castroviejo
Corneoscleral Scissors
small blades
K4-1000 right
K4-1010 left



Castroviejo
Corneoscleral Scissors
medium blades
K4-1020 right
K4-1030 left



Scissors



Troutman-Castroviejo
Corneoscleral Scissors
with stop, small blades
K4-1200 right
K4-1210 left



Troutman-Castroviejo
Corneoscleral Scissors
with stop, medium blades
K4-1220 right
K4-1230 left



Scissors



Vannas

Capsulotomy Scissors
sharp tips, 5mm long blades
angled forward

K4-5020



Vannas

Capsulotomy Scissors
sharp tips, 5mm long blades
angled to side

K4-5030



Gills-Vannas

Capsulotomy Scissors
sharp tips
7mm long blades
straight

K4-5050



Gills-Vannas

Capsulotomy Scissors
sharp tips
7mm long blades
curved

K4-5060



Scissors



Vitreotomy Scissors

20 gauge shank
Stainless steel
construction with
serrated plastic handle.
Autoclavable.



Scissors
horizontal opening
straight
K4-5200



Scissors
horizontal opening
angled 30°
K4-5210



Scissors
vertical opening
angled 60°
K4-5240



Scissors
vertical opening
angled 80°
K4-5250



Katena "Squeeze Handle" Scissors

These "squeeze handle"
instruments feature titanium
handles to reduce overall
weight. 20 gauge shanks
and front ends are made of
hardened stainless steel for
superior sharpness.



Scissors
horizontal, straight
K4-5300



Scissors
horizontal, curved
K4-5310



Scissors
horizontal, angled 45°
K4-5320



Scissors
vertical, 40°
K4-5330



Scissors
vertical, 80°
K4-5340

Eye Scissors
pointed tips
small
K4-7200 straight
K4-7210 curved



Eye Scissors
pointed tips
standard
K4-7400 straight
K4-7410 curved



Eye Scissors
large rings
pointed tips
standard
K4-7440 straight
K4-7450 curved



Quality without Compromise

Scissors



Stitch Scissors
ribbon type
needle points
curved
K4-8000



Stitch Scissors
needle points
light curve
K4-8100



Stevens

Tenotomy Scissors
blunt tips
standard
K4-8500 straight
K4-8510 curved



Stevens

Tenotomy Scissors
ribbon type
blunt tips
K4-8520 straight
K4-8530 curved



Knapp

Strabismus Scissors
blunt tips
K4-9000 straight
K4-9010 curved



Quality without Compromise

Scissors




Strabismus Scissors

ribbon type
blunt tips

K4-9020 straight

K4-9030 curved



Metzenbaum

Dissecting Scissors
gently curved blades
blunt tips

5 3/4" (14.5cm) long

K4-9040



Eye Instruments

Utility Scissors

one pointed and one
blunt tip
5 3/4" (14.5cm) long
straight
K4-9060



Enucleation Scissors

blunt tips
K4-9910 medium curve
K4-9920 strong curve



Quality without Compromise

Botvin

Iris Forceps

For peripheral iridectomies.

K5-1040 serrated tips

K5-1050 1x2 teeth




Gills-Colibri

Forceps

very fine pointed tips
with tying platform

K5-1450



Colibri

Forceps

very delicate
1x2 teeth, 0.12mm

K5-1500



Harms-Colibri

Forceps

very delicate
1x2 teeth, 0.12mm
with tying platform

K5-1520 stainless steel

KT5-1520 titanium



Forceps



Dyson-Colibri

Forceps
very delicate short tips
extremely thin shanks
1x2 teeth, 0.12mm
with tying platform
K5-1522



Colibri

Forceps
standard
1x2 teeth, 0.4mm
with tying platform
K5-1550



Polack

Corneal Suturing Forceps
very delicate
1mm spread
1x2 teeth, 0.12mm
K5-1560




Hofmann-Polack

Corneal Suturing Forceps
with guide pin
1mm spread
1x1 teeth, 0.2mm
K5-1566




Pierse

Corneal Forceps
very delicate
1x1 teeth, 0.1mm
straight
K5-1800



Fechtner

Conjunctiva Forceps
very delicate
ring shaped tips
with tying platform
K5-1820



Castroviejo-Colibri

Forceps
very delicate
1x2 teeth, 0.12mm
K5-1900



Harms-Colibri

Forceps
very delicate
1x2 teeth, 0.12mm
with tying platform
K5-1902




Forceps




Hyde

Corneal Forceps
2x2 teeth, 0.12mm
with tying platform
K5-1920




Bonn

Forceps
very delicate
1x2 teeth, 0.12mm
K5-2000



Bonn

Forceps
1x2 teeth, 0.12mm
with tying platform
K5-2020 stainless steel
KT5-2020 titanium



Thorpe

Corneal Forceps
2x3 teeth, 0.12mm
with tying platform
K5-2100



Dressing

Forceps
serrated, delicate
curved
K5-4010



Tissue

Forceps
1x2 teeth
delicate, straight
K5-4100



Tissue

Forceps
1x2 teeth
delicate, curved
K5-4110



Tissue

Forceps
1x2 teeth
standard, straight
K5-4140



Forceps



Adson

Dressing Forceps
serrated jaws
straight, 5" (125mm) long
delicate
K5-4150



Adson

Dressing Forceps
serrated jaws
straight, 5" (125mm) long
standard
K5-4152



Adson

Tissue Forceps
1x2 teeth
straight, 5" (125mm) long
delicate
K5-4170



Adson

Tissue Forceps
1x2 teeth
straight, 5" (125mm) long
standard
K5-4172



Bayonet
Dressing Forceps
serrated jaws
5 3/4" (146mm) long
K5-4180



Lester
Fixation Forceps
straight
1x2 teeth
K5-4300



Lester
Fixation Forceps
straight
2x3 teeth
K5-4320



Bracken
Fixation Forceps
1x2 teeth
K5-4500



Forceps




Elschnig
Fixation Forceps
1x2 teeth
K5-4550



Troutman
Superior Rectus
Forceps
1x2 teeth
angled
K5-4700



McPherson
Tying Forceps
5mm long smooth jaws
straight
For 8-0 to 11-0 sutures.
K5-5000 stainless steel
KT5-5000 titanium



McPherson
Tying Forceps
5mm long smooth jaws
angled
For 8-0 to 11-0 sutures.
K5-5010 stainless steel
KT5-5010 titanium



Eye Instruments

Kelman-McPherson

Forceps
7.5mm long smooth jaws
angled


K5-5030 stainless steel
KT5-5030 titanium



Haberle-McPherson

Lens Holding Forceps
7.5mm long jaws
with groove, angled
For grasping IOL haptic.


K5-5040



Bechert

Forceps
10mm long smooth jaws
angled

K5-5049



Gaskin

Fragment Forceps
11mm long smooth jaws
angled

K5-5050




Forceps




Kurstin

Flap Stretching Forceps
disc shaped
serrated jaws
*Designed for stretching
the flap to eliminate striae.*
K5-5056



Buratto

LASIK Flap Forceps
disc shaped
serrated jaws
K5-5058



Mendez

Multi-Purpose
LASIK Forceps
spatulated jaws and tips
smooth grasping surfaces
vaulted
*For lifting, dissecting and
grasping the corneal flap.*
K5-5060



Inamura

Capsulorrhexis Forceps
delicate tapered shanks
to seal the incision during
continuous curvilinear
capsulorrhexis
K5-5070




Tennant

Tying Forceps
extra delicate smooth jaws
with 6mm long platform
round handle with guide pin
curved

For 9-0 to 11-0 sutures.

K5-5230 stainless steel

KT5-5230 titanium



Maumenee

Tying Forceps
concave-convex jaws
round handle with guide pin
For 8-0 to 11-0 sutures.

K5-5234 straight

K5-5235 angled



JAWS CROSS SECTION



Harms

Tying Forceps
standard, smooth jaws
For 7-0 to 9-0 sutures.


K5-5400 straight



Harms

Tying Forceps
standard, smooth jaws
For 7-0 to 9-0 sutures.

K5-5410 curved



Forceps




Helveston

Tying Forceps
standard, smooth jaws
For 4-0 to 6-0 sutures.
K5-5420



Fine

Tying Forceps
standard, smooth jaws
For 7-0 to 9-0 sutures.
K5-5550



Barraquer

Cilia Forceps
lightly angled
Textured jaws provide
superior grasping capabilities.
K5-6000



Katena

Cilia Forceps
slanted jaws
gold plated
K5-6020



Eye Instruments

Jewelers

Forceps
delicate, straight
K5-6530 #3



Jewelers

Forceps
delicate, straight
K5-6540 #4



Jewelers

Forceps
extra delicate, straight
K5-6550 #5




Osher

Iris Forceps
very fine textured jaws
angled
K5-6558







Vitreoretinal Forceps
stainless steel with
serrated plastic handle
autoclavable




Vitreoretinal Forceps
smooth jaws
20 gauge, straight
K5-7500



Vitreoretinal Forceps
serrated jaws
20 gauge, straight
K5-7510



Vitreoretinal Forceps
smooth jaws, 1x2 teeth
20 gauge, straight
K5-7520




Vitreoretinal Forceps
smooth jaws, end gripping
20 gauge, straight
K5-7530



Vitreoretinal Forceps
smooth jaws
20 gauge, vertical
K5-7540

Alio
MICS Capsulorrhexis Forceps
23 gauge curved shaft with
micro grasping tips
K5-7651



Quality without Compromise

Forceps



Shepard
Lens Holding Forceps
duckbill shaped jaws
gently curved
K5-8600



Clayman
Lens Holding Forceps
delicate angled jaws
gently curved
K5-8700



Jameson
Muscle Forceps
with slide lock
4 teeth, child size
K5-9000 right
K5-9010 left



Jameson
Muscle Forceps
with slide lock
6 teeth, adult size
K5-9020 right
K5-9030 left



Eye Instruments

Berke
Ptosis Forceps
with slide lock
longitudinally grooved jaws
K5-9300 20mm long
K5-9310 27mm long



Erhardt
Lid Forceps
K5-9400



Downes
Lid Clamp
*Designed for protecting the
eye and holding the lid margin
during upper lid laser
blepharoplasty.*
K5-9410



Lambert
Chalazion Forceps
round, 8mm I.D.
K5-9450



Wies
Chalazion Forceps
grooved upper jaw
K5-9750



Serrefine
serrated jaws
small
1½" (38mm) long
K5-9850



Hartmann
Mosquito Forceps
extra delicate
3¾" (95mm) long
K5-9900 straight
K5-9910 curved



Quality without Compromise

Forceps



Halsted
Mosquito Forceps
delicate
5 ¾" (145mm) long
K5-9940 straight
K5-9950 curved



Baby Jones
Towel Clamp
cross-action
2 ¼" (55mm) long
K5-9960



Backhaus
Towel Clamp
K5-9964 3 ½" (90mm) long
K5-9965 5 ¼" (135mm) long



Eye Instruments

Needle Holders



Castroviejo
Needle Holder
delicate jaws
straight
K6-1500 with lock
K6-1510 without lock




Castroviejo
Needle Holder
delicate jaws
curved
K6-1520 with lock
K6-1530 without lock



Castroviejo
Needle Holder
standard jaws
straight, with lock
K6-1560



Castroviejo
Needle Holder
standard jaws
curved, with lock
K6-1580



Kalt
Needle Holder
delicate smooth jaws
straight, with lock
K6-1800



Drews
Needle Holder
extra delicate jaws
K6-2990 straight, with lock
K6-2991 straight, without lock
K6-2992 curved, with lock
K6-2993 curved, without lock



Barraquer
Needle Holder
short model
very delicate jaws
straight
K6-3000 with lock
K6-3010 without lock



Needle Holders



Barraquer

Needle Holder
short model
very delicate jaws
curved
K6-3020 with lock
K6-3030 without lock



Troutman

Needle Holder
delicate jaws
strongly curved
K6-3250 without lock
K6-3260 with lock



Barraquer

Needle Holder
very delicate jaws
curved
K6-3300 with lock
K6-3310 without lock



KT6-3300 titanium, with lock
KT6-3310 titanium, without lock

Barraquer

Needle Holder
delicate jaws
straight
K6-3500 with lock
K6-3510 without lock



Barraquer

Needle Holder
delicate jaws
curved

K6-3520 with lock
K6-3530 without lock



Barraquer

Needle Holder
standard jaws
curved

K6-3620 with lock
K6-3630 without lock



Anis

Needle Holder
extra delicate jaws
straight

K6-3750 with lock
K6-3751 without lock



Anis

Needle Holder
extra delicate jaws
curved

K6-3752 with lock
K6-3753 without lock



Kerrison

Bone Nibbling Rongeur
9mm opening

K7-1800 size 0, 3.0mm wide

K7-1801 size 1, 3.5mm wide

K7-1802 size 2, 4.0mm wide

K7-1803 size 3, 4.5mm wide



$\frac{2}{3}$

Citelli

Bone Nibbling Rongeur
7mm opening
2.0mm wide

K7-1805



$\frac{2}{3}$

Lacrimal Instruments



Beyer
Rongeur
single action
K7-1810



2/3

Nasal
Speculum
K7-1820 infant
K7-1821 child
K7-1822 adult



2/3

Freer
Periosteal Elevator
one end blunt
other end semi-sharp
4.5mm wide blades
K7-1830



2/3

Eye Instruments



Lacrimal Instruments

Joseph
Periosteal Elevator
gently curved tip
K7-1835



3/4

West
Lacrimal Chisel
straight
K7-1840



3/4

Freer
Lacrimal Chisel
curved
K7-1845



3/4

Quality without Compromise

Lacrimal Instruments

Mallet
lead filled
stainless steel head
K7-1860



$\frac{2}{8}$

Harms
Trabeculotomy Probe
K7-1000 right
K7-1010 left



Infant
Lacrimal Dilator
very fine blunt tip
K7-2000



Muldoon
Lacrimal Dilator
K7-2100



Eye Instruments

VOLGO CARE Irrigating Manipulators

Lu

Irrigating Toric
Lens Manipulator
smooth bulbous shaped tip
end opening
20 gauge thin wall
K7-5780



Alio

MICS Irrigating Fingernail
fingernail shaped tip
1mm diameter port
18 gauge
For Micro Incision Cataract Surgery.
K7-5860



Neuhann

Irrigating Nucleus Manipulator
Y-shaped tip
dual 0.6mm sideports
19 gauge
K7-5864



Rosenwasser

Irrigating Endothelial Stripper
thin, semi-circular tip
end opening irrigation port
20 gauge thin wall
K7-5897



Quality without Compromise


Irrigating Nucleus Choppers



Fukasaku
Small Pupil Snapper Hook
wedge shaped inferior edge
with blunt tip and front notch
for engaging iris
20 gauge thin wall
K7-5873



Rosen
Phaco Splitter
wedge shaped inferior
edge with blunt tip
20 gauge thin wall
K7-5875




Nichamin
Nucleus Chopper
with wedge shaped tip
offset 60°
20 gauge thin wall
K7-5877



Chang
Phaco Chopper
with modified MicroFinger
for horizontal chopping
20 gauge thin wall
K7-5879




Nagahara
Phaco Chopper
with wedge shaped inferior
edge and round blunt tip
20 gauge thin wall
K7-5881



Nagahara
Phaco Chopper
with wedge shaped inferior
edge and round blunt tip
offset 60°
20 gauge thin wall
K7-5883



Giannetti
Nucleus Chopper
and Manipulator
with semi-sharp sides
and concave inferior surface
20 gauge thin wall
K7-5885




Lindstrom
Irrigator
end opening
20 gauge thin wall
K7-5889

Quality without Compromise

Sterilizing Case, medium
autoclavable plastic
with silicone mat
base, lid and mat are perforated
for steam penetration
outside dimensions: 7 $\frac{3}{4}$ " x 4 $\frac{1}{4}$ " x 1"
(200 x 110 x 25mm)
*Accommodates 4 to 6 microsurgical
instruments.*
K9-2025




Sterilizing Case, large
autoclavable plastic
with silicone mat
base, lid and mat are perforated
for steam penetration
outside dimensions:
10 $\frac{3}{8}$ " x 6 $\frac{3}{8}$ " x 1"
(265 x 160 x 25mm)
*Accommodates 10 to 15
microsurgical instruments.*
K9-2030



Sterilizing Cases



Sterilizing Case, double autoclavable plastic double level includes removable tray and two silicone mats base, lid, tray and mats are perforated for steam penetration outside dimensions: 10 $\frac{3}{8}$ " x 6 $\frac{3}{8}$ " x 1 $\frac{7}{8}$ " (265 x 160 x 50mm) Accommodates 20 to 25 microsurgical instruments.
K9-2040




Katena Sterilizing Case, small stainless steel base, lid and mat are perforated for steam penetration outside dimensions: 8 $\frac{1}{4}$ " x 5" x $\frac{3}{4}$ " (210 x 125 x 20mm) Accommodates 8 to 12 microsurgical instruments.
K9-2100 complete
K9-2120 silicone mat only




Eye Instruments


Katena Sterilizing Case, medium stainless steel base, lid and mats are perforated for steam penetration includes 2 silicone mats, rack for 4 cannulas and rack for 4 delicate hooks or knives
outside dimensions: 12 1/2" x 8 3/4" x 3/4"
(320 x 220 x 20mm)
Accommodates 20 to 25 microsurgical instruments.
K9-2200



Katena Sterilizing Case, large stainless steel base, lid and mats are perforated for steam penetration includes 3 silicone mats, racks for 8 cannulas and a rack for 4 delicate hooks or knives
outside dimensions: 14 1/4" x 11 1/2" x 3/4"
(360 x 295 x 20mm)
Accommodates 30 to 40 microsurgical instruments.
K9-2300



Super Blade
sterile, disposable (pk/12)
K20-1760 15°
K20-1770 30°
K20-1780 45°



Chuck Handle
round, serrated
K20-1910 10cm
K20-1920 7.5cm



Bard-Parker Handle
for disposable slot blades
K20-1950 #3
K20-1970 #9



Bard-Parker Blade
sterile, disposable (pk/50)
K20-1981 #11
K20-1985 #15



Performance Microsurgical Blades



K20-1910 handle, 10cm
K20-1920 handle, 7.5cm

Miniature Blades

K20-2464 #64



K20-2465 #65



K20-2467 #67



K20-2469 #69



Lamellar Blade

K20-2466 #66



Sclerotome Blade

K20-2457 #57



Phaco Blade

K20-2475 3.0mm, str.



K20-2476 3.0mm, ang.



K20-2478 3.2mm, ang.



Implant Blade

K20-2488 5.2mm



K20-2490 6.2mm



Super Sharp

K20-2440 15°, 3mm




K20-2442 15°, 5mm



Performance blades are supplied sterile and packaged 5 per box.

Quality without Compromise



Veterinary



Surgical



Beauty Care

VOLGO CARE INTERNATIONAL

Head Office:

Unit No. 4443, DMCC Business Center
Level No.1, Jewellery & Gemplex 3,
P.O. Box 337614 , Dubai, U.A.E.
Tel: +97144572404
Fax: +97144572403
email: info.ae@volgocare.com

Pakistan Office:

Pacca Garah, Kashmir Road, Sialkot Pakistan
Email: info.pk@volgocare.com

UK Office:

Unit 10 Endeavour House ,2 Cambridge Road
Kingston Upon Thames London , KT1 3JU
email: info.uk@volgocare.com

Spain Office:

Badalona, Barcelona Spain
Email: info.sp@volgocare.com

CERTIFIED:

ISO 9001:2015 / ISO13485:2016 